



handicare
Making everyday life easier

OneWaySlide

2075 / 2070

For positioning in a chair or a wheelchair



OneWaySlide is a sliding mat with many functions. The main one improves sitting positions in chairs or wheelchairs. The sliding mat effectively retards forward movement yet slides backwards, making it easy for the carer to help the user move into a more comfortable position in the chair.

OneWaySlide is also available in a tube-shaped design without handles for the more active user.

OneWaySlide is made of polyester, polyethylene and ULF (Ultra Low Friction), Handicare's own unique sliding material.

A non-slide area prevents the user from sliding forward.

Robust handles mean that carers need not physically pull users. Straps for attachment to wheelchairs.

OneWaySlide is available in two models, with and without handles, and in several sizes.



2070 OneWaySlide, tubular 50 x 45 cm



2075 OneWaySlide 130 x 40 cm

For all order enquiries please
email admin@medix21.co.nz or
Tel 04 233 1627

Exclusive to:
Medix 21
29 Jarden Mile
Ngauranga 6035
Wellington

