



Variable pedal resistanceThe adjustable brake can provide drag or lock the wheels for safe transfer.



Adjustable pedals

Backrest anglePosition the backrest vertically or at a 7½° or 15° recline.



drive belt

Carbon fiber drive beltThe drive belt will not rust or stretch, and will last at least five times longer than a bicycle chain.



Seat depthThe seat depth adjusts easily with a knob underneath.

Proudly distributed by:

Medix 21 27 Jarden Mile, Ngauranga, Wellington 6035 New Zealand www.medix21.co.nz admin@medix21.co.nz



