

Support Station rifton®

The Support Station brings new dignity, convenience and possibilities to hygiene care for people with disabilities. In addition to hygiene, it can also be used to aid any sit-to-stand transfer.

With the Support Station, even partially weight-bearing individuals can receive hygiene care in an upright position. This makes the process more dignified, improving access to clothing and skin and eliminating heavy lifting on the part of the caregiver.

Greater independence and dignity: Putting an end to infant-like horizontal hygiene can have dramatic emotional benefits. Even people with limited weight-bearing ability can participate in their own personal care.

Improved safety for caregivers: Eliminating transfers to a mat or changing table means less lifting and struggling with clothing. Better overall accessibility makes the process healthier for everyone.

Practice in strength and skill-building: Using their own strength and skills to transfer and support themselves, people improve weight-



DESIGNED FOR USE WITH THE
MOVE
PROGRAM

Optional Accessories

Support Strap




Kneeboard



Handholds



Support Station dimensions

User dimensions (cm)		K710
Height		117–203
 <p>Important: User height and weight must not exceed the maximum.</p>		
Item dimensions (cm)		K710
Trunk board size		61 x 61
Maximum trunk board height		114
Minimum trunk board height		66
Height of frame		166
Width of frame		70
Depth when trunk board horizontal		79
Depth when collapsed against wall		20
Max. working load (kg)		113

2019 order form

Use dimension chart to select appropriate size.



		K710
* Frame		K715 <input type="checkbox"/>
Support strap		K716 <input type="checkbox"/>
Kneeboard		K717 <input type="checkbox"/>
Handholds	Pair	K718 <input type="checkbox"/>
	Single	K719 <input type="checkbox"/>
Total		K710

* Required

For all order enquiries please
email admin@medix21.co.nz or
Tel 04 233 1627

Exclusive to:
Medix21
29 Jarden Mile
Ngauranga 6035
Wellington