SB6c Roll-InBuddy Lite For roll in showers

The SB6c mobile shower chair from Shower-buddy is designed for accessible or roll-in showers. It can roll straight over most toilets and into the shower, reducing manual handling and transfers for caregivers. The SB6c has a 450mm (18") wide PU foam seat cushion and can also be used as a commode chair.

The SB6c shower chair comes with 4 locking high quality caster wheels, a commode system and numerous accessories to choose from.

The SB6c Roll-InBuddy Lite is lightweight, durable and has a lifetime warranty on the frame



Product features:

- Weight capacity 136kg (300lb)
- Rolls over most standard toilets
- Low maintenance
- Arms removable / lockable / flip-back
- Height adjustable
- Footrests removable and swing-away, height adjustable
- Corrosion resistant Aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts lap or chest belts plus calf strap
- Caster Wheels 100mm (4") all locking (Robby brand)
- Commode system with deflector splash guard
- Tool-free assembly
- Full details see page 31-36

Product Dimensions:

- Length including footrests 900mm (35.4")
- Length without Footrests 700mm (27.5")
- Width over arms 555mm (21.8")
- Width Inside Arms 460mm (18.1")
- Width Between Rear legs 406mm (16")
- Back Width 440mm (17.3")
- Backrest Height from Cushion 350mm (13.7")
- Toilet Height Clearance (maximum) 480mm (18.9")
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 457mm W x 429mm D (18" x 16.9")
- Seat Cushion Height from Floor 500mm 575mm (Adjustable) (19.6" 22.6" adjustable)
- Seat Cushion Opening 180mm W x 330mm L (7" x 13")
- Seat Cushion to Footrest 410mm 510mm (adjustable) (16.1" 20.1" adjustable)









For all order enquiries please email: admin@medix21.co.nz or ph: (04) 233 1627 www.medix21.co.nz

Medix 21 7 Jarden Mile Ngauranga 6035 Wellington



QUICK REFERENCE GUIDE





























4	1
V	

	SB1	SB2	SB2T	SB3T	SB6c	SB6c-22	SB6c-26	SB6w	SB7c	Н5	P1	P2
FEATURES	Shower	Bath Series	Series			Shower	ver			Hybrid	Paediatric	atric
Height Adjustable	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
Lowest height (mm)	550	560	630	510	500	500	500	570	510	500	Note 1	Note 1
Lowest height (inches)	21.6"	21.9"	24.8"	20"	19.6"	19.6"	19.6"	22.4"	20"	19.6"	Note 1	Note 1
Highest height (mm)	620	645	715	600	575	575	575	595	510	575	Note 1	Note 1
Highest height (inches)	24.4"	25.4"	28"	23.6"	22.6"	22.6"	22.6"	23.4"	20"	22.6"	Note 1	Note 1
Commode & Toilet Chair	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Tilt N Space	No	No	Yes	Yes	No	No	No	No	No	No	Note 2	Note 2
Double Locking Struts on Tilt	N.A.	N.A.	Yes	Yes	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	Note 2	Note 2
Transfers into Raised Shower	Yes	2	2	No	No	No	No	No	No	Note 7	Yes	Yes
Transfers into level Shower	2	ł	3	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Note 3	Note 3
Rotates in Shower	Yes	N _O	No	No	Note 4	Note 4	Note 4	Note 4	Note 4	Note 4	Note 5	Note 5
Transfers over Bath	No	Yes	Yes	No	No	No	No	No	N _o	Note 7	Note 6	Note 6
Locking Castor wheels x 4	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	2 Locking	No	N.A.	N.A.
Neck Rest	Accessory	Accessory	Yes	Yes	No	No	No	No	No	No	Yes	Yes
Double safety Brakes	Yes	Yes	Yes	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	Yes	N.A.	N.A.
Lap / Chest Belts	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes	Yes
Soft Seat Cushion	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Soft Overlay / Day Chair Cushion	Accessory	Accessory	Accessory	Accessory	Accessory	No	No	Accessory	No	No	No	No
Foldaway Arms	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
Removable Arms	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
Foldaway Footrest	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes
Removable Footrest	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes	Yes
Tool less easy assembly	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
Non Corrosion Construction	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Weight Capacity Kg	160kg	160kg	160kg	160kg	136kg	170kg	192kg	136kg	120kg	160kg	30kg	90kg
Weight Capacity Lb	350lb	350lb	350lb	350lb	300lb	375lb	425lb	300lb	265lb	350lb	66lb	200lb
2 Year Limited Warranty	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes

but would not normally be applied to	denotes that the chair could be used	
but would not normally be applied to	denotes that the chair could be used	

P1 and P2 accessory is fitted to SB1, SB2, SB2T or SB3T the situation.

Note 1

Note 5 Note 6 Note 7

Note 4

Note 3

so the height is dependant upon which model it is fitted to with the seat height matching the model spec + 25mm (1") P1 and P2 can be tilting or non-tilting depeding upon which model chair the accessory is fitted to. It tilts when fitted to SB2T and SB3T only.

SB3T chair. P1 and P2 can be used to transfer into a level shower when using

Could rotate by turning whole chair provided the shower size allowed this.

P1 and P2 rotate within the shower when fitted to our SB1 chair only.

P1 and P2 will transfer over a bath rub when used in conjunction with SB2 or SB2T chairs

H5 Hybrid is designed for those with some mobilty, yet requiring medium support.

It is anticipated that users could lft the H5 over a shower step or bath tub to use the H5 within the shower or bath tub for support while standing or sitting.