

2026 order form for small

Use dimension chart to select appropriate size.

		Standard R820 small	Hi/lo R830 small
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R821 <input type="checkbox"/> R831 <input type="checkbox"/>
		Adjustable w/ back spring	R822 <input type="checkbox"/> R832 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R823 <input type="checkbox"/> R833 <input type="checkbox"/>
		With spring	R824 <input type="checkbox"/> N/A
* Arm supports (pair)	CHOOSE	Armrests	R827 <input type="checkbox"/> R827 <input type="checkbox"/>
		Forearm prompts	R809 <input type="checkbox"/> R809 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R803 <input type="checkbox"/> N/A
		Long	R894 <input type="checkbox"/> N/A
		Short with casters	R805 <input type="checkbox"/> N/A
		Long with casters	R896 <input type="checkbox"/> N/A
* Pads	CHOOSE	Tan	R815 <input type="checkbox"/> R815 <input type="checkbox"/>
		Blue	R816 <input type="checkbox"/> R816 <input type="checkbox"/>
		Red	R817 <input type="checkbox"/> R817 <input type="checkbox"/>
		Pink	R971 <input type="checkbox"/> R971 <input type="checkbox"/>
		Purple	R972 <input type="checkbox"/> R972 <input type="checkbox"/>
		Green	R973 <input type="checkbox"/> R973 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897 <input type="checkbox"/> R897 <input type="checkbox"/>
		Flat	R898 <input type="checkbox"/> R898 <input type="checkbox"/>
		Adjustable winged	R881 <input type="checkbox"/> R881 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828 <input type="checkbox"/> R828 <input type="checkbox"/>
		Small with chest strap	R819 <input type="checkbox"/> R819 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)		R878 <input type="checkbox"/>	Included
Tray (requires armrests)		R667 <input type="checkbox"/>	R667 <input type="checkbox"/>
Handhold		R649 <input type="checkbox"/>	R649 <input type="checkbox"/>
Chest strap (wide)		R802 <input type="checkbox"/>	R802 <input type="checkbox"/>
Butterfly harness		R834 <input type="checkbox"/>	R834 <input type="checkbox"/>
Thigh belt		R811 <input type="checkbox"/>	R811 <input type="checkbox"/>
Pelvic harness		R845 <input type="checkbox"/>	R845 <input type="checkbox"/>
Hip guides (pair)		R814 <input type="checkbox"/>	R814 <input type="checkbox"/>
Abductor		R884 <input type="checkbox"/>	R884 <input type="checkbox"/>
Adductors (pair)		R892 <input type="checkbox"/>	R892 <input type="checkbox"/>
Leg prompt		R854 <input type="checkbox"/>	R854 <input type="checkbox"/>
Footboard (recommended for chairs with long legs)		R847 <input type="checkbox"/>	Included
Mini kit (backrest insert and trunk support)		R812 <input type="checkbox"/>	R812 <input type="checkbox"/>
Ankle straps (pair)		R899 <input type="checkbox"/>	R899 <input type="checkbox"/>
Sandals (pair)		R944 <input type="checkbox"/>	R944 <input type="checkbox"/>
Wedges (pair, require sandals)		R697 <input type="checkbox"/>	R697 <input type="checkbox"/>
Lumbar and seat support kit		R886 <input type="checkbox"/>	R886 <input type="checkbox"/>
Backrest filler pad		R849 <input type="checkbox"/>	R849 <input type="checkbox"/>
Whitmyer headrest adapter		R804 <input type="checkbox"/>	R804 <input type="checkbox"/>
Universal device mount		R967 <input type="checkbox"/>	R967 <input type="checkbox"/>

* Required

Room to grow

Small chair with mini kit at its smallest dimensions

- 1 Backrest height 32 cm
- 2 Trunk support width 14 cm
- 3 Seat depth 18 cm
- 4 Seat width 18 cm
- 5 Armrest height above seat 13 cm
- 6 Seat height above footboard 18 cm



Small chair at its largest dimensions

- 1 Backrest height 39 cm
- 2 Distance between laterals 29 cm
- 3 Seat depth 30 cm
- 4 Seat width 30 cm
- 5 Armrest height above seat 19 cm
- 6 Seat height above footboard 30 cm



Small standard base
R820 Rifton Activity Chair

User dimensions (cm)	
Height with mini kit	81–122 71–97
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (cm)	R820 small Standard base
Frame width	short legs: 53
	long legs: 58
	short legs w/ casters: 57
	long legs w/ casters: 60
Seat height above floor	short legs: 24–32
	long legs: 47–55
	short legs w/ casters: 34–42
	long legs w/ casters: 47–55
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45° – 110°
Seat height above footboard	18–30
Seat width with hip guides without hip guides	18–23
	30
Seat depth with mini kit	20–30 18–28
Armrest height above seat	13–19
Trunk support width	14–29
Backrest height	32–39
Headrest height above seat	37–53
Max. working load (kg)	34



Small hi/lo base
R830 Rifton Activity Chair

User dimensions (cm)	
Height with mini kit	81–122 71–97
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (cm)	R830 small Hi/lo base
Frame width	66
Seat height above floor	25–64
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	18–30
Seat width with hip guides without hip guides	18–23 30
Seat depth with mini kit	20–30 18–28
Armrest height above seat	13–19
Trunk support width	14–29
Backrest height	32–39
Headrest height above seat	37–53
Max. working load (kg)	34