

2025 order form

Use dimension chart to select appropriate size.

		Z110 small	Z120 medium	Z130 large
* Seat and back (includes seat belt)		Z111 <input type="checkbox"/>	Z121 <input type="checkbox"/>	Z131 <input type="checkbox"/>
* Base	Mobile (tilt-in-space with footboard)	Z112 <input type="checkbox"/>	Z122 <input type="checkbox"/>	Z132 <input type="checkbox"/>
	Mobile (non-tilt with footboard)	Z113 <input type="checkbox"/>	Z123 <input type="checkbox"/>	Z133 <input type="checkbox"/>
	Stationary (non-tilt with footboard)	Z184 <input type="checkbox"/>	Z185 <input type="checkbox"/>	Z186 <input type="checkbox"/>
	Tub base (non-tilt)	Z138 <input type="checkbox"/>	Z138 <input type="checkbox"/>	N/A
	Mounting bar (round-front toilet)	Z115 <input type="checkbox"/>	Z115 <input type="checkbox"/>	Z135 <input type="checkbox"/>
	Mounting bar (elongated toilet)	Z125 <input type="checkbox"/>	Z125 <input type="checkbox"/>	Z135 <input type="checkbox"/>
* Armrests (pair)	Low	Z101 <input type="checkbox"/>	N/A	N/A
	Medium	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>
	High	N/A	Z103 <input type="checkbox"/>	Z103 <input type="checkbox"/>
Casters for mobile base	Small casters (8 cm)	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>
	Large casters (10 cm)	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>
Seat and back pads (recommended)	Rear-open seat pad	Z153 <input type="checkbox"/>	Z163 <input type="checkbox"/>	Z173 <input type="checkbox"/>
	Rear-closed seat pad	Z154 <input type="checkbox"/>	Z164 <input type="checkbox"/>	Z174 <input type="checkbox"/>
	No-hole seat pad	Z187 <input type="checkbox"/>	Z188 <input type="checkbox"/>	Z189 <input type="checkbox"/>
	Soft seat pad	Z924 <input type="checkbox"/>	Z925 <input type="checkbox"/>	Z926 <input type="checkbox"/>
Headrest		Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>
Push handles (pair)		Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>
Butterfly harness		Z116 <input type="checkbox"/>	Z126 <input type="checkbox"/>	Z136 <input type="checkbox"/>
Anterior support/tray		Z145 <input type="checkbox"/>	Z155 <input type="checkbox"/>	Z165 <input type="checkbox"/>
Laterals (pair)	Small	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>
	Large	N/A	Z106 <input type="checkbox"/>	Z106 <input type="checkbox"/>
	Small with strap	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>
	Large with strap	N/A	Z108 <input type="checkbox"/>	Z108 <input type="checkbox"/>
Hip guides (pair)	Small	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>
	Large	N/A	Z149 <input type="checkbox"/>	Z149 <input type="checkbox"/>
	Large with thick pads	N/A	Z198 <input type="checkbox"/>	Z198 <input type="checkbox"/>
Pan		Z117 <input type="checkbox"/>	Z127 <input type="checkbox"/>	Z137 <input type="checkbox"/>
Bowl adapter		Z196 <input type="checkbox"/>	Z196 <input type="checkbox"/>	Z197 <input type="checkbox"/>
Deflector	Standard (requires seat pad)	Z118 <input type="checkbox"/>	Z118 <input type="checkbox"/>	Z128 <input type="checkbox"/>
	Splash guard	Z119 <input type="checkbox"/>	Z119 <input type="checkbox"/>	Z129 <input type="checkbox"/>
	Dynamic deflector	N/A	Z923 <input type="checkbox"/>	Z923 <input type="checkbox"/>
Abductor	Small	Z124 <input type="checkbox"/>	Z124 <input type="checkbox"/>	N/A
	Large	N/A	Z144 <input type="checkbox"/>	Z144 <input type="checkbox"/>
New! Backrest filler pad (Compatible with backrests purchased after February 2026)		N/A	Z922 <input type="checkbox"/>	N/A
Ankle straps (pair, require footboard)		Z156 <input type="checkbox"/>	Z156 <input type="checkbox"/>	Z157 <input type="checkbox"/>
Calf rest (requires footboard)		Z158 <input type="checkbox"/>	Z158 <input type="checkbox"/>	Z159 <input type="checkbox"/>
Footboard (for mounting bar attachment)		Z146 <input type="checkbox"/>	Z147 <input type="checkbox"/>	Z148 <input type="checkbox"/>
Portability base with carry bag		Z139 <input type="checkbox"/>	Z139 <input type="checkbox"/>	N/A


* Required



HTS dimensions



User dimensions (cm)	Z110 small	Z120 medium	Z130 large
Height	76–117	107–142	127–188
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate size of HTS. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>			
Item Dimensions (cm)	Z110 small	Z120 medium	Z130 large
Tilt-in-space mobile base footprint	53 x 69	53 x 69	60 x 84
Tilt-in-space base length at full recline (incl. headrest)	109	125	145.5
Non-tilt base footprint	51 x 64	51 x 64	60 x 74
Seat (front) to footboard height*	22–34	27–39	28.5–44
Seat (front) to floor; mobile tilt frame (add 3 cm for Z181 large casters)	39–55	39–55	44–58
Seat (front) to floor; mobile non-tilt frame (add 3 cm for Z181 large casters)	41–51	43–52	46–56
Seat (front) to floor; stationary non-tilt frame	33–43	36–44	38–48
Seat depth (3 cm increments)	20–28	29–37	37–44
Seat width (depending on hip guides)	20/25/30	20/25/30/36	28/33/38/43
Hole size in seat pad	open 10 x 23, closed 8 x 20	open 11 x 30, closed 9 x 25	open 15 x 36, closed 11 x 30
Armrest height	13.5	17	20
Backrest height above seat	37–44.5	40–19.5	49.5–62
Backrest width	29	34	41
Distance between laterals	12–30	17.5–37	27–46
Tilt in space FORWARD (degrees)	**5°	**5°	**5°
Tilt in space BACK (degrees)	**25°	**25°	**25°
Backrest angle adjustment (degrees)	6° fwd, 6° or 18° back	6° fwd, 6° or 18° back	6° fwd, 6° or 18° back
Footboard angle adjustment (degrees)	10°, 20° or 30° fwd, 10° back	10°, 20° or 30° fwd, 10° back	10°, 20° or 30° fwd, 10° back
Max. clear height under seat (tilt frame)	52	52	52.5
Max. clear height under seat (stationary non-tilt frame)	38	38	42
Max. clear height under seat (mobile non-tilt frame)	46	46	50
Portability base clearance (fits over ADA compliant toilets)	47	47	N/A
Bottom of frame to floor; tilt-in-space mobile base	5	5	5
Tub base footprint	58 x 42 at feet/48 at top	58 x 42 at feet/48 at top	N/A
Tub base height (floor to seat pad)	28 at front, 24 at back	28 at front, 24 at back	N/A
Max. working load (kg)	34	68	***113



*When the mobile base is in the lower positions and the footboard is fully extended, it may prevent the front casters from swiveling. In this case, either change the footboard angle or raise the seat height to give the casters room to swivel.

To achieve 15° of forward tilt for ease of transfer, order the **tilt angle converter, Z169. (The range of adjustment will be 15° back and 15° forward. See reference charts at left.)

***Max. working load of the large footboard is 68 kg.