

Extra wide seat with open sides, 38 cm width for small and medium, 44 cm for the large.

Large size accommodates adults up to 188 cm and 113 kg.

Completely configurable; all sizes of the Rifton Wave fit on any of the 3 base options.

Optional calf rest adjusts independently of seat angle.

Tub transfer base enables users - after being placed in the bath chair using a mechanical lift - to be rolled back and then rotated over the tub for showering.

Conversion kit allows user to economically move up from small to medium size by just purchasing a new back, seat and cover.

Least costly alternative to remodeling a bath room or assisted living, without resorting to bed baths.



Optional Accessories

Head Blocks



Leg Straps



Calf Rest



Height Adjustable Tub Transfer base Legs



Tub Stand



Shower Stand



Tub Transfer Base



Washable Fabric Cover

Pink or Blue



2019 order form

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Select** the color, appropriate bath chair and chest strap.
- 3. Choose** the accessories you need.

The tool-free **conversion kit** converts your small chair to medium as your child grows.

Conversion kit without calf rest

Blue Z268
Pink Z269

(Includes seat, back, cover and chest strap)

Conversion kit with calf rest

Blue Z266
Pink Z267

(Includes seat, back, cover and chest strap)



			Z210 small	Z220 medium	Z230 large
Required components	Chair		Z211 <input type="checkbox"/>	Z221 <input type="checkbox"/>	Z231 <input type="checkbox"/>
	Fabric cover	CHOOSE Blue	Z244 <input type="checkbox"/>	Z246 <input type="checkbox"/>	Z248 <input type="checkbox"/>
		Pink	Z245 <input type="checkbox"/>	Z247 <input type="checkbox"/>	Z249 <input type="checkbox"/>
Chest strap	CHOOSE With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>	
	Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>	
Optional accessories	Head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>
	Leg straps		Z215 <input type="checkbox"/>	Z225 <input type="checkbox"/>	Z235 <input type="checkbox"/>
	Calf rest (calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 <input type="checkbox"/>	Z226 <input type="checkbox"/>	Z236 <input type="checkbox"/>
	Tub stand		Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>
	Shower stand		Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>
	Tub transfer base (calf rest recommended)	CHOOSE Counterclockwise into tub	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>
		Clockwise into tub	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>
	Legs for tub transfer base	CHOOSE Standard	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>
		Long (tub edge over 51 cm high)	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>
	Additional head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>
	Additional fabric cover		Z219 <input type="checkbox"/>	Z229 <input type="checkbox"/>	Z239 <input type="checkbox"/>
Additional chest strap	CHOOSE With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>	
	Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>	
Total		Z210	Z220	Z230	

For all order enquiries please
email admin@medix21.co.nz or
Tel 04 233 1627

Exclusive to:
Medix21
29 Jarden Mile
Ngauranga 6035
Wellington

