

ZING PORTABLE STANDER

PORTABLE , LIGHTWEIGHT, FOLDABLE



Zing

Portable



Supine
Vertical
Prone

The Zing Portable is lightweight (about 10 kg) and features 3 fold-able legs that function like a traditional tripod.

The ergonomic carry handle balances the weight of the unit and makes moving the stander from one place to another easy.

Zing Portable was designed to fold up and pack away for either storage or travel with a minimal footprint.

Portable, lightweight and foldable

- 3 foldable legs, a carrying handle and weight of under 20lbs, makes moving the stander from one room to another a breeze.
- With a small footprint when folded, traveling with your Zing stander is now possible and when a standing session is completed, the Zing Portable can be packed away without taking up much space.
- Great for home, multi-user facility or school.

Plenty of growth

- Fits children from infancy to 36" (91cm), and up to 16kg
- No minimum height or weight requirements which makes starting a standing program as young as possible.
- Comfortable transfer height of 26" (66cm) with near flat loading makes transfers more manageable as children grow.

Easy adjustment

- Built in straps on the support pads.
- Tool-free adjustment of all components.
- Every adjustment can be made with a simple turn of a knob.

Multi use tray supine

- The Supine features a multi-use tray that swings away either direction.
- Tray can remain attached during transfers.
- More activities are possible by swinging tray away during standing therapy.



True hip pivot leg abduction


- Pivot points are in line with the hips, promoting proper alignment in the hips, knees and feet.
- Feet maintain the angle with the rest of the legs - no wedges needed and no additional adjustment to the knees or foot plates.
- Up to 15 degrees of abduction per leg - accomplished through a single knob.
- Promotes musculoskeletal development for children during critical growth years.
- Abduction can be adjusted and used throughout the entire positioning range.

Connect with us



admin@medix21.co.nz | www.medix21.co.nz

Innovative. Portable. Minimal footprint to store



SUPINE

VERTICAL

PRONE



Key user dimensions

Height range	up to 91cm
Weight limit	up to 16.3kg
Weight of basic unit (without tray or mast top)	6.4kg
Weight of Supine Mast Top w/tray and headrest	2.5kg
Weight of Prone Mast Top w/tray	1kg
Tripod Frame Footprint, approximately	91 x 81cm
Transfer Height (Horizontal at hip pivot)	66cm
Leg Abduction	0 to 15°
Foot Plate Size	15 x 6cm
Foot Height Range (Center of Pelvic Support to Footplate)	22-44cm
Kneepad Inside Height, Width is variable (form-to-fit)	10cm
Kneepads Lateral Adjustment, per side	3cm
Kneepads Height Range (Center of pelvic pad to center of knee)	11-34cm
Foot straps- length over arch w/min 1" Velcro overlap	21cm
Tray Size Molded w/8cm Elbow Extensions	19 x 32 cm
Pelvic, Upper Body, and Head Support Pad Sizes	90x 20cm



Connect with us



admin@medix21.co.nz | www.medix21.co.nz