

2026 order form

Use dimension chart to select appropriate size.

Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

Kit without calf rest

- Blue Z268
- Pink Z269

(Kit includes seat, back, cover and chest strap)



Kit with calf rest

- Blue Z266
- Pink Z267

(Kit includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large	
* Chair		Z211 <input type="checkbox"/>	Z221 <input type="checkbox"/>	Z231 <input type="checkbox"/>	
* Fabric cover	CHOOSE	Blue	Z244 <input type="checkbox"/>	Z246 <input type="checkbox"/>	Z248 <input type="checkbox"/>
		Pink	Z245 <input type="checkbox"/>	Z247 <input type="checkbox"/>	Z249 <input type="checkbox"/>
* Chest strap	CHOOSE	With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>
		Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>
Head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>	
Leg straps		Z215 <input type="checkbox"/>	Z225 <input type="checkbox"/>	Z235 <input type="checkbox"/>	
Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 <input type="checkbox"/>	Z226 <input type="checkbox"/>	Z236 <input type="checkbox"/>	
Tub stand		Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>	
Shower stand		Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>	
Tub transfer base (legs required, calf rest recommended)	CHOOSE	Counterclockwise into tub	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>
		Clockwise into tub	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>
Legs for tub transfer base	CHOOSE	Standard	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>
		Long (tub edge over 51 cm high)	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>
Additional head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>	
Additional fabric cover		Z219 <input type="checkbox"/>	Z229 <input type="checkbox"/>	Z239 <input type="checkbox"/>	
Additional chest strap	CHOOSE	With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>
		Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>

* Required

Wave dimensions

User dimensions (cm)	Z210 small	Z220 medium	Z230 large
Height	up to 117	102–142	127–188
 <p>Key user dimension: height Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.</p> <p>Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.</p>			
	Item dimensions (cm)	Z210 small	Z220 medium
Seat depth	27	33	43
Seat width	38	38	44
Back height	58	67	81
Seat height	5–18	5–18	5–18
Seat height (incl. tub stand)	32–47	32–47	32–47
Seat height (incl. shower stand)	70	70	70
Seat height (incl. tub transfer base) with long leg inserts	53–66 62–75	53–66 62–75	53–66 62–75
Length of mesh area when flat (without calf rest)	89	104	128
Length of mesh area when flat (with calf rest)	108	124	155
Base width	46	46	53
Base length (legs down) with legs flat	61 75	61 75	61 75
Tub stand height (legs flat)	20	20	20
Tub stand length (legs flat)	99	99	99
Tub stand width (legs flat)	44	44	44
Shower stand height	56	56	56
Shower stand width	69	69	69
Shower stand length	86	86	86
Tub transfer base footprint	71 x 76	71 x 76	71 x 76
Max. tub height for transfer base with long leg inserts	51 61	51 61	51 61
Maximum working load (kg)	34	68	113
Back and seat angle adjustments (degrees)	0°, 15°, 30°, 45°, 60°, 75°		



All three sizes are available in blue and pink.